

Sleep/Wake Scoring — Installation & Guide

A semi-automated GUI for classifying sleep/wake states from EEG + EMG (and optional video/motion) recordings. A Random Forest model proposes a state for every 4-second epoch; a human reviewer inspects the spectrogram, EMG, theta/delta ratio, and (optionally) video, then corrects the predictions. Corrected labels can be fed back to retrain the model.

This document covers: how to install the package, how to prepare and lay out your data, the configuration file, the full scoring workflow, the GUI controls, and a reference for every file type and data structure the engine consumes or produces.

Package source: `neuroscience_sleep_scoring`

1. What the pipeline does

RAW DATA (per experiment)	PREPROCESSING (run ONCE per experiment)	SCORING (run per acquisition)
AD0_*.mat (EEG ch0) AD2_*.mat (EEG ch2) AD3_*.mat (EMG) EEG/EMG trigger_times.mat vectors *timestamp*.csv (video) prediction *motion*.csv (video) *.mp4 / *.avi	<code>extract_data.py</code> • pick acquisitions • bandpass + downsample • compute normVal • (optional) combine video • (optional) make EDF • (optional) velocity array	<code>New_SWS.py</code> • load downsampled • build feature • Random Forest • review/correct • save StatesAcq*.npy • (optional) retrain

Two scripts you run directly:

Script	When	Frequency
<code>extract_data.py</code>	After acquiring data, before scoring	Once per experiment
<code>New_SWS.py</code>	To score / correct an acquisition	Every scoring session

2. Installation

2.1 Prerequisites

- **Python 3.6+** (3.8–3.10 recommended; see version pitfalls in §2.4).
- **A real GUI display.** The interface uses an interactive Matplotlib backend (`plt.ion()`, `waitforbuttonpress`) plus OpenCV windows for video. It will **not** run on a headless server. Use a local Mac/Windows/Linux desktop, or X-forwarding / a virtual desktop (VNC). On macOS the `MacOSX` or `TkAgg` backend works; on Linux install `python3-tk`.
- **Access to your raw data.** Paths are configured per-experiment (see §4), but a few legacy helper functions still contain hardcoded `/Volumes/...` paths (see §2.4).

2.2 Python dependencies

Declared in `setup.py`:

```
numpy  scipy  matplotlib  psutil  pandas  natsort  opencv-python  pyedflib
scikit-learn  seaborn
```

2.3 Install steps

```
# 1. Clone

git clone https://github.com/YaoChenLabWashU/neuroscience_sleep_scoring.git
cd neuroscience_sleep_scoring

# 2. (Recommended) create an environment

conda create -n sleepscore python=3.9
conda activate sleepscore

# 3. Install the package + dependencies

pip install -e .

# 4. Make the PKA_Sleep dependency importable (see warning below)
```

2.4 ⚠ Paths

Hardcoded NAS paths in a few functions. `New_SWS.build_model`, `SWS_utils.get_eeg_file`, and parts of `extract_data` contain literal `/Volumes/yaochen/Active/...` paths. The main `extract_data` → `New_SWS` flow is driven entirely by the `Score_Settings` file and does **not** depend on these, but the legacy/bootstrap helpers do. Edit them as needed.

3. Required inputs: file types & data structures

3.1 Raw electrophysiology: `AD<channel>_<acq>.mat`

One MATLAB `.mat` file per channel per acquisition, named by channel and acquisition number:

File	Meaning
AD0_<n>.mat	EEG channel 0 (used for the spectrogram and the model)
AD2_<n>.mat	EEG channel 2 (second EEG; enables the 2-channel model)
AD3_<n>.mat	EMG

<n> is the **acquisition number** (an integer, e.g. AD0_4.mat).

Internal structure. Each file contains a single MATLAB struct variable whose name equals the filename stem (e.g. file AD0_4.mat → variable AD0_4). The raw 1-D sample vector is reached through nested indexing, and acquisition metadata lives in a `UserData.headerString` field:

```
mat = scipy.io.loadmat('AD0_4.mat')
data = mat['AD0_4'][0][0][0][0] # 1-D numpy array of samples
hdr = mat['AD0_4']['UserData'].item()['headerString'].item()
fs = int(hdr[hdr.find('inputRate=')+10:].split('\r')[0]) # sampling rate
```

- Sampling rate (`fs`, e.g. 400 Hz) is declared in your config and must match the rig. The pipeline downsamples to `fsd` (e.g. 200 Hz).
- All AD* files for one experiment live in a single directory (`rawdat_dir`).

3.2 Acquisition start times: `trigger_times.mat` (optional but recommended)

If present in `rawdat_dir`, this file provides exact start datetimes for each acquisition (variable `trigger_times`, indexed by acquisition order). It is used to align video/motion to ephys.

- **If absent**, the code falls back to estimating each acquisition's start from the `.mat` file's filesystem modification time minus its duration. This is less reliable; prefer providing `trigger_times.mat` whenever video is used.

3.3 Video: <basename><n>.mp4 (or .avi)

Optional. Required only if `vid = 1`. One video per recording segment, named <basename><n>.<ext> (e.g. 1tFLiPAKARmemEEGEMG0095_3.mp4). Stored in `video_dir`. `.mp4` is searched first, then `.avi`.

3.4 Bonsai timestamp CSVs: `*timestamp<n>.csv`

Optional (needed for video and/or motion). One CSV per video segment, stored in `csv_dir`. Format: a single column of ISO-8601 frame timestamps, no header:

```
2023-04-12T15:04:01.2345678+00:00
2023-04-12T15:04:01.2680000+00:00
...
```

The <n> suffix in the filename ties each timestamp file to its video and motion file. Timestamps are auto-adjusted against the file's modification time to correct clock offset.

3.5 Motion CSVs: `*motion<n>.csv`

Optional (needed only if `movement = 1`). One per segment, in `csv_dir`. Two supported formats, selected by the `DLC` config flag:

- **DeepLabCut (`DLC = 1`)** — standard multi-row DLC output. The body part named by `DLC Label` (e.g. `center`) is read using its `_x`, `_y`, and `_likelihood` columns. Frames with `likelihood < 0.8` are flagged as low-quality.
- **Bonsai raw (`DLC = 0`)** — a headerless CSV with either `x, y` columns or `Timestamps, X, Y` columns.

Velocity per epoch is computed as the Euclidean displacement of the tracked point between epoch boundaries.

3.6 Directory layout (recommended)

```
<rawdat_dir>/                                # all raw ephys for one experiment
  AD0_4.mat  AD0_5.mat  ...
  AD2_4.mat  AD2_5.mat  ...
  AD3_4.mat  AD3_5.mat  ...
  trigger_times.mat                                # optional

<video_dir>/                                  # if vid = 1
  <basename>3.mp4  <basename>4.mp4  ...

<csv_dir>/                                    # if vid/movement = 1
  *timestamp3.csv  *timestamp4.csv  ...
  *motion3.csv     *motion4.csv     ...

<savedir>/                                    # created by extract_data.py
  (downsampled data + scoring outputs — see §6)

<model_dir>/                                  # the shared model + logs
  EEG_2chan_EMG_nomovement.joblib
  mouse_2chan_model.pkl
  *_scoringlog.txt
```

4. Configuration: `Score_Settings.json`

Everything is driven by a single JSON file passed as the only command-line argument. Keep a personal copy per experiment. Every key:

Key	Type	Meaning
<code>basename</code>	string	Experiment name; used to build <code>video/normVal</code> filenames.
<code>rawdat_dir</code>	path	Folder containing the raw <code>AD*.mat</code> (and <code>trigger_times.mat</code>).

Key	Type	Meaning
model_dir	path	Folder holding the shared .joblib model and *_model.pkl.
video_dir	path	Folder with video files (if vid).
csv_dir	path	Folder with Bonsai timestamp/motion CSVs.
modellog_dir	path	Where the model scoring log is written (usually = model_dir).
personallog_dir	path	Where your personal scoring log (personal_scoringlog.csv) is written.
species	string	"mouse"
mod_name	string	Model name stem (e.g. "mouse_2chan"). Determines *_model.pkl filename.
mouse_name	string	Animal ID (not the experiment name).
epochlen	int	Epoch/bin length in seconds. Minimum 4 ; the model and features assume 4.
fs	int	Raw sampling rate of the .mat data (e.g. 400).
fsd	int	Target downsampled rate (e.g. 200).
emg	0/1	1 if EMG (AD3) was recorded.
vid	0/1	1 to enable video review.
movement	0/1	1 to use motion/velocity as a feature & plot.
EEG channel	list[int]	EEG channels to use, e.g. [0, 2]. Must be a list. Length 2 → 2-channel model.
EMG channel	int	The EMG channel number (e.g. 3).
Acquisition	list[int]	Acquisitions to process. Auto-filled by extract_data.py.
Filter High	int	Band-pass high cutoff (Hz) for EEG/EMG (e.g. 100).
Filter Low	float	Band-pass low cutoff (Hz) (e.g. 0.5).
savendir	path	Where downsampled data and scoring outputs go.
Bonsai Version	int	Which Bonsai workflow produced the video/CSV files.
DLC	0/1	1 if motion CSVs are DeepLabCut output; 0 for raw Bonsai.
DLC Label	string	Body part to track for velocity (e.g. "center").
Minimum_Frequency	int	Lower bound (Hz) of the displayed spectrogram (e.g. 1).
Maximum_Frequency	int	Upper bound (Hz) of the displayed spectrogram (e.g. 20).
vmin	number or "None"	Spectrogram color floor (dB). "None" = auto.
vmax	number or "None"	Spectrogram color ceiling (dB). "None" = auto.

5. Usage workflow

Step 1 — Preprocess

```
cd neuroscience_sleep_scoring
python extract_data.py /path/to/Score_Settings.json
```

This runs, in order:

1. **choosing_acquisition** — lists the `ADO_*.mat` files and asks whether to add **all** acquisitions, choose them **individually**, or include only those above a **minimum length**. Your selection is written back into the config's `Acquisition` list.
2. **downsample_filter** — band-pass filters (Butterworth, order 3, using `Filter Low-Filter High`) and resamples each acquisition from `fs` to `fsd`. EMG is additionally 10 Hz high-passed. Saves the downsampled arrays (see §6.1).
3. **get_normalizing_value** — computes the median total spectral power across acquisitions per EEG channel and saves it as `<basename>_normVal.npy`. This value normalizes band-power features so they are comparable across sessions.
4. **EDF export** (prompted, optional) — bundles EEG+EMG into `.edf` files for external viewers.
5. **Video/motion** (only if `movement = 1`) — `combine_bonsai_data` merges all timestamp and motion CSVs into `All_timestamps.pkl / All_movement.pkl`, and optionally builds the full `velocity_vector.npy`.

Step 2 — Score / correct (per acquisition)

```
python New_SWS.py /path/to/Score_Settings.json
```

The engine then:

1. Lists available acquisitions and which already have a saved `State` file.
 2. Prompts: **which acquisition** to score.
 3. Splits the acquisition into ≤ 1 -hour segments (most are a single segment).
 4. For each segment, prompts: **check existing scoring (c)** or **score a new dataset (s)**:
 - **s (score new)**: optionally runs the Random Forest to predict states, saves `model_prediction_Acq<a>_hr<h>.npy`, then opens the GUI pre-filled with the prediction for you to correct. (Choosing not to use the model opens a blank, all-zero scoring you fill in manually.)
 - **c (check/fix)**: loads the existing `StatesAcq<a>_hr<h>.npy` and opens the GUI so you can adjust it. If any epoch is still unscored (state 0) it warns you and offers to go back.
 5. On exit, prompts whether to **update the model** (retrain with this segment's corrected labels) and whether to **update your personal log**.
-

6. The GUI: layout & controls

When scoring opens, you get two figure windows:

Figure 1 (the main figure) — 5 stacked panels:

1. Spectrogram of EEG channel 0 (frequency band set by `Minimum/Maximum_Frequency`).
2. **Predicted states** strip (the row you click to correct) with a black model-confidence trace overlaid.
3. Spectrogram of EEG channel 2.
4. Velocity (or "No Movement Available").
5. EMG amplitude.

Figure 2 (the zoomed/raw figure) — synchronized detail panels for the EMG, theta/delta ratio, and (if available) velocity around the current epoch, plus a moving marker line.

Controls

Action	How
Re-center / inspect an epoch	Click once inside the spectrogram (panel 1). The raw-trace figure (and video, if enabled) jumps to that epoch $\pm$ surrounding epochs.
Correct a run of epochs	Click the start epoch in the predicted-states strip (panel 2), then click the end epoch. Switch to the terminal and type the new state number (1–4) when prompted. The strip recolors and <code>StatesAcq*.npy</code> is re-saved immediately.
Toggle a crosshair line	With Figure 2 focused, press <code>l</code> to drop/erase a vertical alignment line across panels.
Score a single epoch (manual mode)	Keys <code>1,2,3,4</code> set the current epoch's state.
Exit a video clip early	Press <code>v</code> .
Finish this segment	Press <code>d</code> (done). You'll then be asked about saving / updating the model.
Quit immediately	Press <code>q</code> (closes everything and exits).

All keys are lowercase (no Shift). If a click "does nothing," make sure Figure 2 (not the terminal or video window) has focus. Don't toggle the line while hovering the motion panel.

7. Sleep-state encoding

States are integers. In the GUI they are entered as keystrokes 1–4 and color-coded in the predicted-states strip:

State	Color	Conventional meaning*
0	grey / white	Unscored (placeholder; must be resolved before retraining)
1	green	REM / active sleep
2	blue	NREM / slow-wave sleep
3	red	Wake
4	purple	Fourth state (e.g. active vs. quiet wake / transitional)

8. Data-structure reference

8.1 Generated by `extract_data.py` (written to `savendir`)

File	Shape / type	Description
<code>downsampEEG_Acq<a>.npy</code>	1-D float array	Full downsampled EEG for acquisition a (used to get total length).
<code>downsampEMG_Acq<a>.npy</code>	1-D float array	Full downsampled EMG for acquisition a.
<code>AD<chan>_downsampled/downsampEEG_Acq<a>_hr<h>.npy</code>	1-D float array	Per-hour EEG segment for channel <code>chan</code> , trimmed to a whole number of epochs.
<code>downsampEMG_Acq<a>_hr<h>.npy</code>	1-D float array	Per-hour EMG segment.
<code>AD<chan>_downsampled/<basename>_normVal.npy</code>	scalar	Median total power; normalizes band-power features.
<code>All_timestamps.pkl</code>	pandas DataFrame	Columns <code>Timestamps</code> , <code>Filename</code> — all video frame times.

File	Shape / type	Description
All_movement.pkl	pandas DataFrame	Columns Timestamps, X, Y, Likelihood, Filename — all tracked positions.
velocity_vector.npy	2×N array	Row 0 = per-bin velocity, row 1 = bin time (s).
*.edf	EDF	Optional EEG+EMG export for external viewers.

8.2 Generated by scoring (New_SWS.py, written to savedir)

File	Shape / type	Description
StatesAcq<a>_hr<h>.npy	1-D int array, one entry per epoch	The primary output — your final/working sleep-state labels. Re-saved every time you correct a run of bins.
model_prediction_Acq<a>_hr<h>.npy	1-D int array	The raw Random Forest prediction (before your corrections), saved when you score in model mode.

8.3 Model

File	Type	Description
<mod_name>_model.pkl	pickled pandas DataFrame	Accumulated training table : every feature column plus State, animal_name. Grows each time you "update the model".

8.4 Logs

File	Type	Description
<mod_name>_scoringlog.txt	text	Who scored what, when, and in which mode (corrected / scored with ML model / scored in legacy mode).
personal_scoringlog.csv	CSV	Columns Date, Mouse Name, Acquisition, State Array Location.

8.5 In-memory structures

d (settings dict) — the parsed `Score_Settings.json`, threaded through nearly every function.

FeatureDict — built by `SWS_utils.build_feature_dict`. One key per feature, each a 1-D array of length `num_epochs`. Keys are generated per EEG channel `c` (e.g. `EEGChannel0`, `EEGChannel2`):

- `EMGvar` — EMG variance per epoch (if EMG present).
 - `<c>var` — EEG variance per epoch.
 - `<c>Delta`, `<c>Theta`, `<c>Alpha`, `<c>BroadTheta`, `<c>Fire` — band powers (0.5–4, 5–8, 8–12, 2–16, 4–20 Hz), normalized by `normVal`.
 - `<c>nb` — narrow-band theta = `Theta / BroadTheta`.
 - `<c>thet_delt` — `Theta / Delta` ratio.
 - `<c>delta_pre/post, ..._pre2/post2, ..._pre3/post3` — the delta feature shifted $\pm 1/2/3$ epochs (gives the model temporal context). Same for `theta_*` and `nb_pre/post`.
 - `Velocity` — per-epoch velocity (if movement).
 - `State` — added after scoring; the supervised target.
-

9. Module / function map

What each file is responsible for, and the key functions the engine calls.

`New_SWS.py` — scoring engine (entry point)

- `load_data_for_sw(config) → start_swscoring(d)` — top-level driver.
- `start_swscoring(d)` — acquisition/segment loop; prompts; orchestrates predict → review → save → optional retrain.
- `display_and_fix_scoring(...)` — builds the figures, wires up the cursor, runs the interactive event loop, writes `StatesAcq*.npy`.
- `update_model(d, FeatureDict)` — appends corrected labels to the training table and retrains the `.joblib`.
- `model_log(...), personal_log(...)` — write the audit logs.

`extract_data.py` — preprocessing (run once)

- `choosing_acquisition, pulling_acqs` — discover/select acquisitions.
- `downsample_filter` — filter + resample → downsampled `.npy`.
- `get_normalizing_value` — compute `normVal`.
- `combine_bonsai_data, make_full_velocity_array` — video/motion consolidation.
- `make_edf_file, save_to_edf, get_phys_fs` — optional EDF export & header parsing.

`SWS_utils.py` — shared library (~50 functions)

- **Features:** `build_feature_dict`, `generate_signal`, `bandPower`, `plot_spectrogram/my_specgram`, `post_pre`, `get_ThD`, `get_total_power`, `prepare_feature_data`, `get_xfeatures`, `adjust_movement`.
- **Model:** `update_sleep_df`, `update_sleep_model`, `retrain_model`, `random_forest_classifier`, `build_joblib_name`, `load_joblib`, `fix_states`, `define_microarousals`.
- **Figures/interaction:** `create_prediction_figure`, `plot_predicted`, `create_zoomed_fig`, `plot_zoomed_data`, `update_raw_trace`, `make_marker`, `add_buffer`, `clear_bins`, `correct_bins`, `print_instructions`.
- **Video/timestamps/motion:** `load_video`, `pull_up_movie`, `get_videofn_from_csv`, `timestamp_extracting`, `pulling_timestamp`, `initialize_vid_and_move`, `movement_extracting`, `movement_processing`, `get_movement_segs`, `sort_files`, `sort_timestamp_files`, `get_AcqStart`.
- **DLC:** `movement_extracting` (DLC branch), `DLC_check_fig`, `rename_DLC_csvs`, `transfer_DLC_files`.

SW_Cursor.py — interaction classes

- `Cursor` — full interaction: re-center clicks on the spectrogram, two-click bin selection on the states strip, crosshair, and the d/l key handlers. Sets flags (`replot`, `change_bins`, `DONE`) that the `New_SWS` event loop reads.
- `ScoringCursor` — a trimmed variant.

10. Troubleshooting

Symptom	Likely cause / fix
<code>ModuleNotFoundError: PKA_Sleep on import</code>	The undeclared dependency isn't on the path. Install/add <code>PKA_Sleep</code>
<code>AttributeError: module 'numpy' has no attribute 'float' after an odd keypress</code>	<code>NumPy ≥ 1.24</code> . Pin <code>numpy < 1.24</code> or patch <code>np.float('nan') → float('nan')</code>
<code>ImportError: cannot import name 'simpson'</code>	SciPy too old. Upgrade SciPy or alias <code>simps</code>
Figure opens but clicks do nothing	Wrong window focused — click into Figure 2 , not the terminal/video.

Symptom

Likely cause / fix

"No videos found" / timestamp errors

Check `video_dir/csv_dir`, filename suffixes (`<basename><n>`, `*timestamp<n>`, `*motion<n>`), and that counts of timestamp vs motion files match.

Wrong/old labels reload when checking

You're loading `StatesAcq<a>_hr<h>.npy` — confirm the acquisition/hour and that prior corrections were saved (they save automatically on each bin edit).

Model not found in score mode

The `.joblib` whose name matches your EEG channel/emg/movement config isn't in `model_dir`.
